


NCAC CNAC

National Companion Animal Coalition - La coalition nationale sur les animaux de compagnie



UPDATED LIST / MISE À JOUR 2026 - 01- 01

A review process has been established by the National Companion Animal Coalition (NCAC) to assess the conformity of radio frequency identification (RFID) products and processes with the revised Canadian standard for electronic (microchip) identification of companion animals. The revised Canadian standard came into effect December 1, 2004 and is based on ISO microchips, accompanying databases and compatible readers to promote maximum long-term recovery potential of companion animals. Only those RFID products that have been submitted for review and have been deemed in conformity with the Canadian standard will be recognized as suitable for companion animal identification. Only ISO microchips are acceptable under the Canadian standard. However, the NCAC requests companies to ensure backwards compatibility of readers and database capability to accept the older standard until the year 2020. In addition, the NCAC permits users to apply non-Canadian standard microchips only if required and specially requested for export purposes.

Un processus d'examen a été établi par la Coalition nationale pour les animaux de compagnie (CNAC) afin d'évaluer la conformité des produits et des processus d'identification par radiofréquence (IRF) par rapport à la norme canadienne révisée pour l'identification électronique (par puce) des animaux de compagnie. La norme canadienne révisée est entrée en vigueur le 1^{er} décembre 2004 et respecte la norme des puces de l'ISO, des bases de données correspondantes et des lecteurs compatibles afin de faciliter le potentiel de récupération à long terme des animaux de compagnie. Seuls les produits IRF qui ont été soumis aux fins d'examen et ont été jugés conformes à la norme canadienne seront reconnus comme étant acceptables pour l'identification des animaux de compagnie. Seules les puces de l'ISO sont acceptables en vertu de la norme canadienne. Cependant, la CNAC demande aux compagnies de garantir la compatibilité rétroactive des lecteurs et la capacité des bases de données d'accepter l'ancienne norme jusqu'en 2020. De plus, la CNAC autorise les utilisateurs à employer des puces standard non canadiennes seulement en cas de besoin et uniquement si elles ont fait l'objet d'une demande spéciale aux fins d'exportation.

NCAC CNAC

National Companion Animal Coalition - La coalition nationale sur les animaux de compagnie



NCAC Statement on Implantation of Multiple Transponders in Animals

Background

The National Companion Animal Coalition (NCAC) is aware of circumstances where an animal may be previously implanted with a microchip (transponder) that is not recognized by the NCAC (see the NCAC recognized which can be found on CVMA's website at:

<http://canadianveterinarians.net>)

Some of these non-recognized microchips can be non-ISO or encrypted and may not be read by some readers, which can result in unsuccessful lost pet recovery.

NCAC Recommendation

In order to help ensure successful lost pet recovery, the NCAC recommends that an NCAC-recognized ISO microchip be implanted. The NCAC does not recommend that the non-recognized microchip be removed and rather recommends the owner contact the non-recognized microchip company to update the owner information in their database. This is so the database information will match that of the new NCAC-recognized microchip that is implanted, and the correct owner information is accessed no matter which microchip is scanned.

Énoncé de la CNAC portant sur l'implantation de multiples micropuces dans un animal

Mise en situation

La Coalition nationale pour les animaux de compagnie (CNAC) est consciente qu'un animal puisse avoir été implanté avec une micro puce (transpondeur) non reconnu par la CNAC (voir la liste de produits reconnus par la CNAC sur le site web de l'ACMV :

<http://canadianveterinarians.net>)

Ces micropuces non reconnues peuvent être de technologie non-ISO ou même cryptée, les rendant illisibles par certains lecteurs, affectant l'efficacité du processus de recouvrement.

Recommandation de la CNAC

Afin d'aider à assurer un recouvrement efficace, la CNAC recommande qu'une micropuce de technologie ISO reconnue par la CNAC soit implantée. La CNAC ne recommande pas que l'autre micro puce soit enlevée mais recommande plutôt que le propriétaire de l'animal contacte la compagnie ayant vendue la micro puce non reconnue afin que celle-ci ajoute cette nouvelle information à leur base de données.

De cette manière, la bonne information sera accessible peu importe la micro puce qui est lue.

NCAC CNAC

National Companion Animal Coalition - La coalition nationale sur les animaux de compagnie



Recognized RFID Products

Reconnaissance des produits à identification par radiofréquence

| Company <i>Fabricant/Distributors</i> | Microchip Model <i>Modèle de puce (Injectable)</i> | Reader Name/Model <i>Nom/modèle du décodeur</i> | Animal Identification Database <i>Base de données d'identification des animaux</i> |
|--|--|--|--|
| BC SPCA 1-855-622-7722 1-866-303-5950 https://bcpetregistry.ca/ | ISO FDX-B Manufacturer Code 933 | TN009B-PRD Multitag TN Pocket Tracker | BC Pet Registry 1-855-622-7722 info@bcpetregistry.ca |
| Canadian Kennel Club (416) 674-3699 1-855-364-7252 www.ckc.ca | Trovan ID 162 ISA FDX-B Manufacturer code: 956 | Petscan Reader | Canadian Kennel Club 1-855-364-7252 1-800-396-1896 www.ckc.ca |
| Datamars Inc. 1-781-281-2216 www.datamars.com www.petlink.net | Datamars TIP 8010 Datamars T-IT8100 Datamars T-IS 8010 Datamars Microfindr Slim Datamars Microfindr Value ISO FDX-B <i>Manufacturer code: 981</i> | CompactMax Universal Scanner ISO Standards 1178/45 | Petlink International database 1-877-738-5465 www.petlink.net |
| EIDAP Inc. 1-888-346-8899 905-970-8181 www.eidap.com | Trovan ID 162 ISO FDX-B <i>Manufacturer code: 956</i> AEG ID 162 ISO FDX-B <i>Manufacturer code: 968</i> | Trovan LID560 <i>Pet-Identifier * Pocket-Reader EX* Portable Scan* Reader RT100V2* * No longer Marketed</i> | EIDAP Inc. 1-888-346-8899 905-970-8181 www.eidap.com |
| Insight Microchips Corporation 1-888-664-9002 www.insightmicrochips.com | Insight Microchips Manufacturer code (900) Manufacturer Code Range: Current Range: 219.000.000.000 to 219.000.999.999 Updated Range: 219.000.000.000 to 219.002.999.999 | Insight Microchips | PrimeTrackr Microchips https://www.insightmicrochips.com/ |
| Intervet Canada Corp 1-888-466-3242 www.HomeAgain.com | HomeAgain® TempScan™ ISO FDC-B Manufacturer code: 985 | Universal WorldScan™ Reader Plus | HomeAgain 1-888-466-3242 www.HomaAgain.com |



NCAC CNAC

National Companion Animal Coalition - La coalition nationale sur les animaux de compagnie




| | | | |
|---|---|---|---|
| <p>Microchip 4 Solutions Inc. 1-877-738-4384 (519) 348-8491 www.microchipsolutions.com</p> | <p>ISO FDX-B <i>Manufacturer code: 952</i></p> <p>ISO FDX-B <i>Manufacturer Code: 939</i></p> | <p>PetScan V2 PetScan V8</p> | <p>Microchip 4 Solutions Inc. 1-877-738-4384 (519) 348-8491 www.microchipsolutions.com</p> |
| <p>Pethealth Inc. 1-866-597-2424 (905) 842-2615 www.24petwatch.com</p> | <p>Allflex Injectable ISO FDX-B</p> <p>Allflex Injectable ISO FDX-B Mini-Chip Manufacturer code: 982</p> <p>Allflex Europe SAS 982061 FDX-B ISO TCHIP-CDN 24 Petwatch. Manufacturer Code: 982</p> <p>24 Petwatch. Manufacturer: Datamars, Manufacturer code (941)</p> | <p>Allflex Portable</p> <p>Allflex ScanFlex AFX-100 Reader</p> <p>Destron pocket Reader EX</p> <p>N/A</p> | <p>Pethealth Inc. 1-866-597-2424 (905) 842-2615 www.24petwatch.com</p> |
| <p>Petidco / Avid Canada The Microchip Company 1-800-338-1397 (403) 264-6300 www.petidco.com</p> | <p>Identichip ISO FDX-B Manufacturer code: 965</p> <p>ISO FDX-B Manufacturer code: 900 (shared) Identification Code range: 900.026.000.000.000 to 900.026.000.999.999</p> | <p>N/A</p> | <p>Petidco / Avid Canada The Microchip Company 1-800-338-1397 (403) 264-6300 www.petidco.com</p> |
| <p>Identrac Inc. 1 877 906 3030 www.identrac.ca</p> | <p>EasyTrac-ID ISO FDX-B Manufacturer code: 967</p> <p>Smartrac ISO FDX-B Manufacturer code: 992</p> | <p>Smartrac ADS Intrace handheld reader (RH5)</p> | <p>Identrac Inc. 1-877-906-3030 www.identrac.ca</p> |
| <p>SyrVet Canada (450) 361-1504 1-888-779-7838</p> | <p>AEG ID ISO FDX-B Manufacturer Code: 968</p> | <p>AEG ID ARE H4 AEG ID ARE H5</p> | <p>aZooome Identification 1-844-429-6663 www.azoo.me</p> |



NCAC CNAC

National Companion Animal Coalition - La coalition nationale sur les animaux de compagnie



| | | | |
|--|--|---|--|
| <p>Viaguard Inc 416-691-4167 1-877-842-4827 www.accu-metrics.com</p> | <p>ISO FDX-B Product Code: 900056 Manufacturer Code: 900 (shared) Identification Code starts 990 000000 7...</p> | <p>N/A</p> | <p>Viaguard Inc 1-877-842-4827 www.accu-metrics.com</p> |
| <p>Microchip ID Systems Inc. (985) 898-0811 1-800-434-2843 www.microchipIDsystems.com</p> | <p>ISO FDX-B Buddy ID Mini-Chip <i>Manufacturer code: 933</i></p> | <p>N/A</p> | <p>Microchip Registration Centre (for pets) 1-800-434-2843 www.BuddyID.com</p> |
| <p>Vetygo +1-866-438-5067 (438) 506-7665 vetygo.ca</p> | <p>ISO FDX-B Manufacturer code 991001 inclusive ISO HT003 FDX-B Manufacturer code 991012 Transponder (microchip) model: ISO FDX-B Manufacturer code 991003 inclusive Transponder (microchip) model: ISO FDX-B Manufacturer code (900) Range (217.000.000.000 217.000.999.999)</p> | <p>Animal ID PT160 Animal ID PT200 N/A N/A</p> | <p>Vetygo By email: info@vetygo.ca 866-438-5067</p> |

